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Guidance and manufacture's declaration - electromagnetic emissions

The device is suitable for use in the specified electromagnetic environment and it has meets the following standard's emission requirements.

Phenomenon	Profession healthcare facility environment	Home healthcare environment
Home healthcare environment	CISPR 11, Group 1, Class A or B	CISPR 11, Group 1, Class B
Harmonic distortion	IEC 61000-3-2, Class A or not applicable	NA
Voltage fluctuations and flicker	IEC 61000-3-3 or not applicable	NA

Guidance and manufacture's declaration - electromagnetic immunity

The device is suitable for use in the specified electromagnetic environment and it has meets the following immunity test levels. Higher immunity levels may cause the device's essential performance lost or degraded.

Phenomenon	Basic EMC standard or test method	Professional healthcare facility environment	Home healthcare facility environment
Electrostatic discharge	IEC 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	
Radiated RF EM fields	IEC 61000-4-3	3V/m 80MHz-2.7GHz 80%AM at 1kHz or 2Hz	10V/m 80MHz-2.7GHz 80%AM at 1kHz or 2Hz
		1kHz or 2Hz can be specified by the manufacturer	

Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See the RF wireless communication equipment table in "Recommended minimum separation distances".
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m; 50 Hz or 60Hz
Electric fast transients bursts	IEC 61000-4-4	NA
	For input a.c. power port d.c. power lines or signal input/output lines whose length exceeding 3m	
Surges	IEC 61000-4-5	NA
	For 1. input a.c. power port; 2. all d.c. power ports connected permanently to cables >3m; 3. output signal output lines connected directly to outdoor cables	
Conducted disturbances induced by RF fields	IEC 61000-4-6	NA
	For 1. input a.c. power port; 2. all d.c. power ports connected permanently to cables >3m; 3. all patient-coupled cables; 4. SIP/SOP whose maximum cable length $\geq 3m$	
Voltage dips	IEC 61000-4-11	NA
Voltage interruptions	IEC 61000-4-11	NA
UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz or 30 cycles at 60Hz		

Recommended minimum separation distances

Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or system's basic safety and essential performance may be affected. This device has been tested with the immunity level test and listed below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this device as recommended below.

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM $\pm$ 5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						

Huai An Hua Lve Electronics Co., Ltd.
No.211 Yingchun Avenue, Xucheng Industrial Zone, Xuyi,
Huai'an City, Jiangsu 211700, China

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg Germany
E-Mail: shholding@hotmail.com
Tel.: 049-40-2513175 Fax.: 049-40-255726

E-Invasion A/S
Artillerivej 86, 1
2300 Copenhagen
Denmark
E-mail: info@babyplan.com